



# DHVI Biomolecular Interaction Analysis (BIA) Core Facility

## Mission


The DHVI Biomolecular Interaction Analysis (BIA) Core Facility offers state-of-the-art **Surface Plasmon Resonance (SPR)** and **Biolayer Interferometry (BLI)** instruments, expertise, training, and support for **label-free, real-time** analysis of biomolecular interactions to basic and clinical researchers within the Duke Community.

## Leadership / Experience


- Director: S. Munir Alam, PhD**
  - Munir Alam is a Professor of Medicine and Pathology at Duke University, a DHVI Principal Investigator, and the Founding Director of the DHVI BIA Core Facility. Dr. Alam is a recognized expert in biochemical and biophysical analyses of receptor-ligand interactions, including TCR-ligand and antibody-antigen binding, and HIV-1 Envelope protein antigenicity characterization. He is Director of Protein Analytics at both the DHVI and its GMP facility, a Co-Investigator in the Duke CHAVD program, and was the Proteomics Director for the Duke CHAVI program.
- Manager: Kenneth Cronin**
  - Ken Cronin has served as Lab Manager for Dr. Alam for 7 years and recently has taken on role of BIA Core Manager. He has extensive experience in assay development and biophysical characterization of molecular interactions by SPR, BLI, and FPLC.

## Services

- Binding Specificity and Screening**
  - Yes/No Binding
  - Ranking of Binding Response
- Kinetics and Affinity**
  - Direct Binding Kinetics
  - Steady State Affinity
  - In-solution Affinity
  - Avidity Assessment
- Active Concentration**
  - Surface Competition
  - Solution Inhibition
- Thermodynamics**
- Epitope Mapping**
- Conformation Changes**




## Technology / Instruments / Resources




- Biacore SPR systems are **low sample consumption** (typically < 25  $\mu$ g) and are **compatible with a wide-range of biomolecules**
  - Proteins, peptides, DNA/RNA, small molecules, carbohydrates, polymers, lipid membranes, micelles, liposomes, polymers, whole viruses, and whole cells
  - Crude cell culture lysates or clinical sera/bodily fluids
- Biomolecules can be functionalized on sensor chips via:
  - Covalent coupling**
  - High-affinity capture** (with any affinity tag)
  - Hydrophobic adsorption**

## DHVI BIA Core Facility Platforms




**Biacore 3000**

- Standard platform
- Simultaneously detects 4 interactions
- 10 RU Sensitivity




**Biacore T200/S200**

- LMW screening
- Simultaneously detects 4 interactions
- 1 RU Sensitivity
- T200 available for Independent Use



**ForteBio Octet-RED96**


- Biolayer Interferometry (BLI)
- Ligand-coated biosensors in 96-well plate
- Crude sample compatibility
- Simultaneously detects 8 interactions



**ForteBio Octet HTX**

- High-throughput Biolayer Interferometry (BLI) analysis Supports 96- and 384-well plate racks

## Recent Projects / Publications



- Kinetics/Affinity** – Henderson R., Watts, B.E., et al., *Selection of immunoglobulin elbow region mutations impacts interdomain conformational flexibility in HIV-1 broadly neutralizing antibodies*. Nature Communications **2019**, 10 (654).; Enam, S.F., et al., *Enrichment of endogenous fractalkine and anti-inflammatory cells via aptamer-functionalized hydrogels*. Biomaterials **2017**, 142.
- Vaccine-induced Avidity Screening** – Lynch, H.E., et al., *Surface plasmon resonance measurements of plasma antibody avidity during primary and secondary responses to anthrax protective antigen*. Journal of Immunological Methods **2014**, 404.
- Competition** – Tomaras, G.D., et al., *Vaccine-induced plasma IgA specific for the C1 region of the HIV-1 envelope blocks binding and effector function of IgG*. PNAS-USA **2013**, 110 (22).
- Epitope Mapping** – Liao, H.-X., et al., *Vaccine Induction of Antibodies against a Structurally Heterogeneous Site of Immune Pressure within HIV-1 Envelope Protein Variable Regions 1 and 2*. Immunity **2013**, 38 (1).
- Nanoparticle/Liposome Antigenicity** – Alam, S.M., et al., *Role of HIV membrane in neutralization by two broadly neutralizing antibodies*. PNAS-USA **2009**, 106 (48).

## Reservations / Service Requests

- Independent Use** – Self-Service, T200 only
  - Complete training in the independent operation of the Biacore T200
  - Available 24/7/365 through CoreResearch@Duke
- Sample Submission** – Limited Full-Service, All platforms
  - All experimentation, optimization, and analysis will be performed by Core Staff.
- Consultation**
  - Assistance in experimental design and analysis/interpretation of data
- Contact Ken Cronin to schedule an initial consultation

## Contact Us

- BIA Core Facility website:**  
[dhvi.duke.edu/programs-and-centers/shared-resources/cores/biomolecular-interaction-analysis-bia-core-facility](http://dhvi.duke.edu/programs-and-centers/shared-resources/cores/biomolecular-interaction-analysis-bia-core-facility)
- Email:** kenneth.cronin@duke.edu
- Phone:** (919) 684-4380