

Manual Gating of Lymphocytes using BD MultiSET Software

Prepared by Sylvester Hood for Immunology Quality Assessment

IQA Program Manager: Raul Louzao

In this presentation, you will learn how to adjust:

- CD45 gating
- Bead gating
- CD3/SSC gating

We will be using BD MultiSET software, but the same principles may be applied to other flow cytometry gating programs and platforms.



MultiSET File Edit Tools Help Mon 12:27 PM

untitled

Sign In Set Up FACSComp Test Prefs Samples Acquisition Analysis Lab Report Phys Report Summary

Message: Please enter the Sample information and then click the Run Tests button. Sample Name, Sample ID, and Panel Name are required for each entry.

	Sample Name	Sample ID	Case Number	Panel Name	WBC Count (x1000)	Lymphs (%)	Lymphs (0
1	UK NEQAS 212202	UK NEQAS 358	PSE 9 PSM 34_36	4 Color TBNK + TruC	0.0	0.0	0.00
2	UK NEQAS 212202	UK NEQAS 359	PSE 9 PSM 34_36	4 Color TBNK + TruC	0.0	0.0	0.00
3					0.0	0.0	0.00
4					0.0	0.0	0.00
5					0.0	0.0	0.00
6					0.0	0.0	0.00
7					0.0	0.0	0.00
8					0.0	0.0	0.00
9					0.0	0.0	0.00
10					0.0	0.0	0.00
11					0.0	0.0	0.00
12					0.0	0.0	0.00
13					0.0	0.0	0.00
14					0.0	0.0	0.00
15					0.0	0.0	0.00
16					0.0	0.0	0.00
17					0.0	0.0	0.00
18					0.0	0.0	0.00
19					0.0	0.0	0.00
20					0.0	0.0	0.00
21					0.0	0.0	0.00
22					0.0	0.0	0.00
23					0.0	0.0	0.00
24					0.0	0.0	0.00
25					0.0	0.0	0.00
26					0.0	0.0	0.00
27					0.0	0.0	0.00
28					0.0	0.0	0.00
29					0.0	0.0	0.00
30					0.0	0.0	0.00
31					0.0	0.0	0.00
32					0.0	0.0	0.00
33					0.0	0.0	0.00
34					0.0	0.0	0.00
35					0.0	0.0	0.00
36					0.0	0.0	0.00
37					0.0	0.0	0.00

1. Click here to add samples

2. Click here to select a panel to use

3. Click here to start manually gating

The Sample list
Note: You will need to use files ending in a .01 extension. Files ending in .lab or .phy are not used here.

Stop Add Samples... Run Tests


CD3/CD8/CD45/CD4 TruC

Data Set [1] Data File: UK NEQAS 35806.01

Reagent Lot ID: 78222 Events Acquired: 15000 Abs Cnt Bd Lot ID: 20071 Attr Def File: 3/8/45/4 MLT/TruC v2.0

File ID: 6B849D7F-52A6-476F-A1C7-AD01208F7719 Beads Per Pellet: 47100

Lymph Events	3000	
Bead Events	2841	
CD3+ %Lymph	85.97	Hi
CD3+ Abs Cnt	855	Lo
CD3+CD8+ %Lymph	35.00	
CD3+CD8+ Abs Cnt	348	
CD3+CD4+ %Lymph	32.87	
CD3+CD4+ Abs Cnt	327	Lo
CD3+CD4+CD8+ %Lymph	0.27	
CD3+CD4+CD8+ Abs Cnt	3	
CD45+ Abs Cnt	995	
TH/S Ratio	0.94	



QC Messages:

- Code 3: Sample integrity questionable. Visually inspect.
- Code 4: The CD3+ %Lymph value lies outside the normal reference range.
- Code 4: The CD3+ Abs Cnt value lies outside the normal reference range.
- Code 4: The CD3+CD4+ Abs Cnt value lies outside the normal reference range.

Tube Name: CD3/CD16+56/CD45/CD19 TruC

Status: Not yet run

File Name:

Select "Manual Gate"

Stop Rerun Manual Gate Skip... Continue






Lymph	
497	219
691	329
750	306
803	176
689	0
567	0
497	48

Mes	

Manual Gate






Manual Gating Work Table:

1. CD45/SSC
2. CD3/CD8
3. CD3/SSC
4. CD8/CD4




- As samples degrade, monocytes may present lower than usual on the CD45/SSC plot. The automated gating may not always detect this anomaly.
- Double click on the CD45/SSC plot to enlarge it. Redraw the lymphocyte gate if necessary.






- You may observe events from the bead gate in the CD8 gate. This is an indication that you need to tighten the bead gate.



Before



After



- As samples degrade, the CD3+ population will appear dimmer. It may even start to overlap with the CD3- population. The automated software may not always detect this anomaly.
- Double click on the CD3/SSC plot to enlarge it. Adjust the CD3 gate if necessary.






Lymph	
498	130
763	262
803	112
687	4
567	0
498	50

Manual Gate

When you are ready to check your work, click the “Analyze” button for an updated report. Then, re-enter the Manual Gating Worktable by selecting “Manual Gate”







Checking the CD45/SSC Plot

The bead population is completely purple. If there is any white, the bead gate on the CD3/CD8 plot is pulled too high.

The entire lymphocyte population is included. If you see a section that is 90% "blue" events at the top of the gate, you may need a tighter gate. The gate pictured here is well-adjusted.



Checking the CD8/CD4 Plot

- There are few or no bead events in the CD8 gate.
- The CD4+ and CD8+ events are in tight separate clusters.



Checking the CD3 Gating using the CD3/CD8 Plot

- The CD4 and CD8 clusters should be tight and separate. CD3- “Blue” events should not be present in the CD3+ space, and CD4 and CD8 events should not be in the CD3- space. A line break is included below to help visualize.



The staff at the IQA Center is available for technical consultation concerning laboratory techniques for these procedures. You may contact Raul Louzao at 919-684-5861 raul.louzao@duke.edu or Sylvester Hood at 919-613-4469 sfh7@duke.edu

